



# Remains Processor

for Cremated Animals



- Built for rapid, efficient processing of cremated remains.
- Processes the remains down to 1/16" to a fine powder.
- Safe - Push button switch operates the processor.
- When released the processor stops.
- Average processing time: 35 pound animal in 5-10 seconds; 150 pound animal 10-15 seconds.
- Processing chamber is 3/8" thick x 10" round x 13-1/2" tall steel plate for long life with minimal maintenance.
- Cam lock lid with gasket for a complete seal & handle for tilting processor for discharge.
- Unique twin flail blade design with easy access for blade replacement. Blades are 1/4" case hardened steel.
- Blade life - 800+ processes. Replacement blades are economical.
- Space requirements: 20" wide x 36" deep x 44" high. Weight of processor is 220 pounds.
- Electrical Service: 1 hp motor @ 115 volts, 20 amps, single phase or 1-1/2 hp motor @ 230 volts, 15 amps, single phase.
- Supplied with 8' heavy duty electrical cord with NEMA L6-15p plug. Matching receptacle L6-15R required (By Others).